
**INSTRUCTIONS FOR USE**  
MULTIPARAMETER CABLE

-  Extreme temperatures for storage and use
-  Protect from humidity
-  Lot Code
-  Quantity
-  Expiry
-  Medical Device
-  Do not put in the trash. Comply with the European directive on electrical equipment (WEEE).
-   Precautions for use. The instructions for use of the multiparameter cable contain important safety information. Consult the package leaflet for instructions for use
-  Manufacturer
-  Manufacture Date
-  No latex  
Not made with natural rubber latex
-  Non sterile
-  Reference
-  CE marking - EU MDR 2017/745

**MULTIPARAMETER CABLE**

**I. IDENTIFICATION / SCOPE**

The Nissha Medical Technologies multiparameter cable is the ideal intermediary between a patient and an electromedical instrument; it makes it easier to connect the various devices (sensors, electrodes, etc.) used for medium- and/or long-term monitoring of a patient's physiological parameters.

Nissha Medical Technologies offers a full line of multiparameter cables for the monitoring of the following physiological parameters:

- electrocardiogram **ECG**
- temperature **T°**
- blood oxygen saturation **SpO<sub>2</sub>**
- non-invasive pressure **NIP**

To re-order, use the code number marked on your Nissha Medical Technologies multiparameter cable or its packaging.

**Contact Nissha Medical Technologies for more information about this product.**

**II. STORAGE / PACKAGING**

Nissha Medical Technologies multiparameter cables are individually packaged.

When not in use, a **cable** must be stored in its original packaging to prevent any deterioration of its protective sheath, which might impair its service life, its performance, or its level of safety.

Storage conditions are as follows:

- **Temperature: -40 to +70 °C**

**III. PERFORMANCE / RELIABILITY / SAFETY / COMPATIBILITY / MECHANICAL INTEGRITY**

**PERFORMANCE / RELIABILITY:**

Nissha Medical Technologies multiparameter cables are inspected during and after production in accordance with a technical protocol established in conformity with the rules and directives currently in force.  
(/G-MED technical report no. 600 12053-5056 44)

**SAFETY:**

Nissha Medical Technologies multiparameter cables are designed and manufactured in conformity with the specifications of the applicable European and/or national standards currently in force:

**IEC 60601-1/IEC-60601-...(specific to each parameter)**

**SPECIFICATIONS:**

Level of protection from electric shocks:

**see SV  
Class II**

Classification according to the recommendations of EU MDR 2017/745 class I The low-frequency leakage currents measured in accordance with the recommendations of current standards applicable to this product are smaller than the allowed values.

(G-MED technical report no. 600 12053-5056 44)

**COMPATIBILITY:**

To ensure compatibility between the devices:

- Use Nissha Medical Technologies multiparameter cables only with "Nissha Medical Technologies " accessories (ECG conductors, electrodes, sensors, etc.). (see §IV)
- Connect only the accessories and/or devices necessary for installing the Nissha Medical Technologies multiparameter cable mentioned in the commercial documentation of the product concerned. (see §IV)

**MECHANICAL INTEGRITY:**

To enable Nissha Medical Technologies multiparameter cables (conductors, connectors, cable, etc.) to withstand tensile and bending stresses and reduce the risk of deterioration in service, Nissha Medical Technologies uses highly reliable, high-quality materials.

The encapsulated connectors have flexible sleeves to minimize the risk of failure of the cable there.

Nissha Medical Technologies multiparameter cables are designed to withstand repeated defibrillation shocks. They have no accessible metallic parts.

**IV. INSTALLATION / USE / MAINTENANCE / HYGIENE****INSTALLATION / USE:**


For optimal installation and use of the Nissha Medical Technologies multiparameter cable, proceed as follows:  
(see also the directions for use of the electromedical device and of the accessories concerned)

Patient end:



- Place the various sensors and/or electrodes on the chosen patient sites as appropriate for the measurement.
- Connect the intermediate cables from the various sensors and/or conductors (ECG leads, etc.) to the interface (patient end) of the Nissha Medical Technologies multiparameter cable

**Attention:**

Refer also to the directions for use specific to each sensor (SpO<sub>2</sub>, temperature, NIP, etc.) and conductor (ECG leads) connected to the Nissha Medical Technologies multiparameter cable.


- Attach the multiparameter cable near the patient using a suitable attachment system on the underside of the cable.

**Instrument end:**

- Connect the plug(s) of the Nissha Medical Technologies multiparameter cable to the appropriate medical device.

- Power up the electromedical device and check that it is working properly (refer to the directions for use of the device).

**PREVENTIVE MAINTENANCE:**

Before each use:

- check that the insulation (protective sheath of connecting cable, cable interface, connectors) has not suffered any mechanical damage.
- check (at least once a month) the electrical continuity of the conductors using Nissha Medical Technologies physiological signal simulator(s)

(contact Nissha Medical Technologies to know which testers are best suited to your needs)

**GENERAL CONDITIONS OF SERVICE:**

The conditions of service of the multiparameter cable are as follows:

- Ambient temperature: +10 to +40 °C

**HYGIENE:**

When the Nissha Medical Technologies multiparameter cable is used in wards where a certain level of sterility is required, it must be cleaned and/or disinfected as described below.

- Clean the Nissha Medical Technologies multiparameter cable and its connecting cable (connector included) using a cloth soaked with soapy water

**Warning:**

Do not let the liquid infiltrate the connectors of the interface cable.

- Disinfecting the Nissha Medical Technologies multiparameter cable:

**- Method A**

(based on ANIOS Laboratories study no. 6416.94/0387)

- fill a dip tank with a 0.5% solution of HEXANIOS G+R
- partially immerse the cable, protecting its ends to prevent any electrical problem with the connectors
- leave for 15 minutes
- rinse the cable
- dry the cable with absorbent paper
- during the immersion, rub the ends of the cable with an ANIOS LINGET.

(change the solution every 48 hours)

**- Method B** (new low-temperature method)

(based on ANIOS Laboratories study no. 14496.02/052)

- fill a dip tank with activated ANIOXYDE 1000 preparation
- wait 30 minutes
- partially immerse the cable, protecting the ends to avoid any electrical problem with the connectors
- leave for from 10 to 30 minutes, depending on desired level of disinfection
- check the peracetic acid content regularly using the checking strip
- rinse with ANIOS network water (pH = 7.3, TH = 48°F)
- dry the cable with absorbent paper

(see ANIOS documentation for further information)

**IV. WARRANTY / LIABILITY**

Nissha Medical Technologies guarantees for one year any unused cable that is kept in its original packaging and has suffered no visible damage.

Nissha Medical Technologies guarantees the conformity of the device to the specifications of the safety and performance standards currently in force applicable to it.

**ATTENTION :**

The class and type of protection from electric shocks (BF, CF, etc.) are determined by those of the electromedical device to which the multiparameter cable is connected.

Consult the directions for use of the device concerned, and of the accessories to be used, before putting the cable into service. Nissha Medical Technologies cannot be held liable for problems that occur in the event of failure to follow the rules for installation and use given in this document.

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[HS.NisshaMedical.com/IFU](https://www.hs-nissha.com/IFU)

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### CAVO MULTIPARAMETRO

	Extreme temperatures for storage and use			Precautions for use. The instructions for use of the multiparameter cable contain important safety information. Consult the package leaflet for instructions for use
	Protect from humidity			
	Lot Code			
	Quantity			Manufacturer
	Expiry			Manufacture Date
	Medical Device			No latex
	Do not put in the trash. Comply with the European directive on electrical equipment (WEEE).			Non sterile
				Reference
				CE marking - EU MDR 2017/745

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[HS.NisshaMedical.com/IFU](http://HS.NisshaMedical.com/IFU)

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-  Extreme temperatures for storage and use
-  Protect from humidity
-  Lot Code
-  Quantity
-  Expiry
-  Medical Device
-  Do not put in the trash. Comply with the European directive on electrical equipment (WEEE).
-  Precautions for use. The instructions for use of the multiparameter cable contain important safety information. Consult the package leaflet for instructions for use
-  Manufacturer
-  Manufacture Date
-  No latex  
Not made with natural rubber latex
-  Non sterile
-  Reference
-  CE marking - EU MDR 2017/745

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