## INSTRUCTIONS FOR USE

### DISPOSABLE TUBING

-  **CAUTION**
-  **TEMPERATURE LIMITS**
-  **KEEP AWAY FROM DIRECT SUNLIGHT**
-  **PROTECT FROM MOISTURE**
-  **LOT CODE**
-  **QUANTITY**
-  **EXPIRY**
-  **DO NOT USE IF PACKAGING IS DAMAGED**
-  **MEDICAL DEVICE**
-  **DISPOSE OF AFTER 24 HOURS OF USE**
-  **NON-STERILE**
-  **MANUFACTURER**
-  **LATEX FREE**  
Not made with natural rubber latex
-  **PHT FREE**
-  **READ INSTRUCTIONS FOR USE**
-  **REFERENCE**
-  **EC ACCORDING TO DIRECTIVE 93/42/EEC AND UPDATE 2007/47/EC**  
0197

-  Using the tubing for more than 24 hours can cause bacteria to develop in the stagnant water within the tube.
- Before using this device, please read the instructions for the endoscope you will be using, as well as for any accessories that will be used.

## DISPOSABLE TUBING

Please read the following information carefully.

Failure to use the product as intended may have harmful consequences for the patient.

### Important reminder:

This document contains instructions for using the "ENDOtube for Wash Bottle" tubing for endoscopes. It contains no references to particular medical procedures. The manufacturer will not be held liable for any problems caused by improper use of the device.

## I. INDICATION

ENDOtube for Wash Bottle is designed to supply an endoscope with air and water during gastrointestinal endoscopy procedures. ENDOtube for Wash Bottle is a tube for 24 hours of use that connects bottle of sterile water to an endoscope with its direct connector.


There are four models:

- **FLO000001** - DISPOSABLE TUBING - OLYMPUS GI ENDOSCOPES
- **FLP000002** - DISPOSABLE TUBING - PENTAX GI ENDOSCOPES
- **FLF000003** - DISPOSABLE TUBING - FUJIFILM GI ENDOSCOPES 5-600 SERIES
- **FLF000004** - DISPOSABLE TUBING - FUJIFILM GI ENDOSCOPES 700 SERIES

**ENDOtube for Wash Bottle+CO2** is designed to supply an endoscope with air or carbon dioxide (CO2) and water during endoscopy procedures. **ENDOtube for Wash Bottle + CO2** is a tube for 24 hours of use that connects bottle of sterile water to an endoscope with its direct connector.

There are four models:

- **FLO001CO2** - DISPOSABLE TUBING + CO2 OPTION - OLYMPUS GI ENDOSCOPES
- **FLP002CO2** - DISPOSABLE TUBING + CO2 OPTION - PENTAX GI ENDOSCOPES
- **FLF003CO2** - DISPOSABLE TUBING + CO2 OPTION - FUJIFILM GI ENDOSCOPES 5-600 SERIES
- **FLF004CO2** - DISPOSABLE TUBING + CO2 OPTION - FUJIFILM GI ENDOSCOPES 700 SERIES

### CONTRAINDICATIONS:

The contraindications are those particular to any endoscopy procedure.

## II. STORAGE / PACKAGING

### STORAGE:

- Room temperature: +10°C to +50°C
- Expiry date: 2 years (see on the label)

### PACKAGING:

Wash bottle tubing for endoscopes is packaged individually in boxes of 20.  
Once removed from its packaging, the tubing must be used within 24 hours of opening.

## III. WARNINGS AND PRECAUTIONS

### The ENDOfube for Wash Bottle range:

- Must be used within 24 hours of opening.
- Must be disposed of as waste in compliance with the country's applicable regulations.
- Must never be connected or reconnected to an endoscope that has not been washed as per manufacturer instructions.
- Not designed to be reused after being used for 24 hours.
- Must only be used by qualified endoscopy personnel.

Keep an eye on the level of water in the bottle during the procedure. The water bottle must be changed when the water level falls below the tube.

**When changing the water bottle, disconnecting the system between two patients, or during any other procedures, authorised personnel must use appropriate techniques to avoid contaminating the ENDOfube for Wash Bottle system.**  
Be sure to identify any leakages in the connector to avoid wetting the floor and causing staff to slip.

The device comes with a **male Luer** lock fitting.  
Special care must be taken to avoid inadvertently connecting the connector to an existing line on the patient with a female connector.

### COMPATIBILITY:

- FLO000001 - OLYMPUS GI Endoscopes
- FLP000002 - PENTAX GI Endoscopes
- FLF000003 - FUJIFILM GI Endoscopes 5-600 SERIES
- FLF000004 - FUJIFILM GI Endoscopes 700 SERIES
  
- FLO001CO2 - OLYMPUS GI Endoscopes
- FLP002CO2 - PENTAX GI Endoscopes
- FLF003CO2 - FUJIFILM GI Endoscopes 5-600 SERIES
- FLF004CO2 - FUJIFILM GI Endoscopes 700 SERIES

## IV. USE / MAINTENANCE / CLEANING / DISINFECTION / STERILISATION

### USE:

It is important that you familiarise yourself with the device and complete a trial run before using it in a clinical context.  
(See also the instructions for the endoscope and any accessories used)

1. Open a bottle of sterile water (25cl, 50cl, or 100cl).
2. Open the packaging and remove the ENDOfube for Wash Bottle system.
3. Insert the tube into the bottle of sterile water and firmly press the stop to make sure it is watertight.
4. Check that the clamp is in the open position.
5. Connect the air/water connector to the endoscope.
6. Turn on the processor's light source.
7. Before starting the procedure, check that it is working properly by starting the endoscope's air/water canal. If there is a

malfunction, check that the water bottle is properly sealed and properly connected to the endoscope.

8. If you need to replace the water bottle, use suitable techniques.
9. Once the procedure is complete, clamp and unplug the ENDOTube for Wash Bottle system.
10. You can then turn off the processor once you have unplugged the ENDOTube for Wash Bottle system.
11. Make sure that the device is properly disposed of within 24 hours of opening the packaging. A traceability label is provided for this purpose.

**MAINTENANCE:**






- Do not use the device if the packaging is damaged.
- No repairs can be carried out on this device.

**CLEANING / DISINFECTION / STERILISATION:**

It is a non-sterile device for 24-hour use that cannot be sterilized. Disinfection or cleaning is not necessary.

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