



LONG-TERM /ICU PATIENT MONITORING

CLARAVUE is an efficient and hygienic solution for simple, accurate and fully radiotransparent patient monitoring. CLARAVUE includes prewired disposable electrodes for use with patented, reusable bioactive trunk cables.

For the past ten years, CLARAVUE has been known to provide very good trace quality and infection prevention in prominent hospitals globally.

Hygienic, Single Patient, Prewired Electrodes for Long-Term Monitoring

FEATURES

- 

LARGE ELECTRODE AND GEL AREA
A larger surface area adheres better to patient's skin, especially in challenging conditions.
- 

FOAM BACKING
Sturdy material retains integrity during long-term wear.
- 


HYGIENIC, DISPOSABLE SOLUTION
Reduces the risk of cross contamination.
- 

PREWIRED ELECTRODES
Reduces motion artifact and increases patient comfort.
- 

PATENTED CONNECTOR CABLE
Standardised for all monitors.
- 

VISIBLE BACKING LINER
Prevents slip hazard.
- 

BIOACTIVE CABLE
Molybdenum Trioxide coating provides antimicrobial effects and prevents cross contamination.



PACKAGING

Electrodes are packaged in a single-use pouch with an easy-to-open tear.

Single-Use Kit

Box




CLARVUE® RADIOTRANSSPARENT DISPOSABLE PREWIRED ELECTRODES


Product #	Description	Per Kit	Kits/Box	Kits/Case		Solid Gel	Duration of Use	Radiotranslucent
32031769	LPM M5 RT SET (SIDE)	1	20	240				
32031770	LPM M5 RT SET (TOP)	1	20	240				

ELECTRODE ANATOMY

FRONT



BACK



COMPATIBLE MONITOR CONNECTIONS

Compatible cables include adaptors for Bard, Biosense, Mindray, Datascope, GE Marquette, Philips, Datex, Mennen, Mortara, Nihon Kohden, GE, Getemed, Spacelabs, Siemens, Drager, Schiller, Zoll and more.



Bard



Philips



GE Prucka



St. Jude




GE Marquette Hellige



Nihon Kohden

CABLE CONNECTIONS



3-13-01-078 (08/2020)