



NICU | PICU

CLARAVUE® provides the best-in-class, hygienic solution for simple, accurate and fully radiolucent patient monitoring. CLARAVUE® includes pre-wired and pre-positioned disposable electrodes and patented reusable bioactive trunk cables.

For the past ten years, CLARAVUE® has performed as the proven choice for superior trace quality and infection prevention in prominent European hospitals.

Lightweight, disposable, pre-wired electrodes for the NICU | PICU

FEATURES



5-DAY WEAR TIME

Fewer disruptions ensure the body focuses on essential growth and repair.



LIGHTWEIGHT, DISPOSABLE LEADWIRES

Use on even the tiniest of patients. Disposable leadwires reduce the risk of cross contamination.




SOFT & GENTLE HYDROGEL

Provides comfortable contact without the use of harsh acrylic adhesives.



PRE-WIRED ELECTRODES

Help to reduce motion artifact.




PACKAGING

Electrodes are packaged in a single-use pouch with an easy-to-open tear

Box




Single-Use Kit



FRONT

BACK





CLARVUE® RADIOLUCENT DISPOSABLE PRE-WIRED ELECTRODES

Product #	Description	Per Kit	Kits/Box	Kits/Case	Solid Gel	Duration of Use
32028776	Pediatric Patient Kit - 3 electrodes	3	40	480		

COMPATIBLE MONITOR CONNECTIONS

Compatible cables include adaptors for Bard, Biosense, Mindray, Datascope, GE Marquette, Philips, Datex, Mennen, Mortara, Nihon Kohden, GE, Getemed, Spacelabs, Siemens, Drager, Schiller, Zoll and more.



Mindray



Philips



Spacelabs TruLink



GE



Physio



Nihon Kohden

The OEM trademarks identified herein are the trademarks of the respective OEMs, and not Vermed. Vermed disclaims any affiliation, connection, or association between its products and those of the respective OEMs.

CABLE CONNECTIONS

ELECTRODE → **PATENTED ELECTRODE CABLE CONNECTOR** → **CABLE** → **MONITOR**

