

Sustainability Meets Efficiency

Designed to reduce environmental impact through reductions in carbon emissions during production and shipping, EcoVue's proprietary FlexPac[®] is lightweight and requires less energy to manufacture and transport. It is ergonomically designed for ease of use and functionality and allows for an unmatched 99.5% product evacuation rate (versus the standard 85.5% evacuation with most bottles), dramatically reducing waste and driving bottom line gains.

- ✔ 99% natural formulation
- ✔ No dyes, fragrances, or parabens
- ✔ Propylene glycol free to reduce wear and tear on probes
- ✔ Exceptional image quality
- ✔ Water-based formula for easy clean up
- ✔ Lower environmental impact



285 - EcoVue[®] Ultrasound Gel 250g FlexPac[®]



Designed specifically for the 250g FlexPac[®]

Acclimate[®]: Safe Gel Warming for Patient Comfort

Custom-designed for our 250g FlexPac[®] flexible packaging to provide consistent warming of the gel throughout the package. The front can also accommodate single-use packets and the notched cutout allows for easy access to insert and remove them.

The Acclimate[®] warmer provides the convenience of single-button control with a built-in two stage thermal cut-off to protect from overheating. Functions as a wall mount or countertop unit, it is UL listed, available in 120V or 230V, and includes a two-year warranty.

Acclimate[®] helps to deliver a better experience for today's ultrasound providers AND their patients.