



PerformancePlus™ Electrodes

Stress Test Kit

FEATURES



V1-V2

Small rectangular electrode with two sensors



V3-V6:


Small teardrop shaped electrode for placement in tight areas



RA, LA, RR, LL:

Large round electrodes for superior adhesion

Unique single use pouch containing 10 electrodes is perfect for diagnostic stress testing.



PACKAGING

HYGIENIC, DISPOSABLE SINGLE-USE PACKS

Reduces risk of cross contamination of unused electrodes



CASE

FDA / GS1 Compliant Labels



Our unique single use pouch containing 10 electrodes is perfect for diagnostic stress testing. Electrodes come in varying sizes for ease of placement in designated positions.

A10015-10S - STRESS TEST KIT - FOR 10 LEAD PLACEMENT
10 ELECTRODES/ KIT | 50 KITS / CASE

Size	Shape	Per Pouch	Wet Gel
2.25" / 57.15mm	Round	4	
1.5" x 38.1mm	Teardrop	4	
1.75" x 1.625" / 44.43mm x 41.275mm	Square, dual element	1	

WHAT'S IN A KIT?

