



# SilveRest™ Electrodes

Adult Diagnostic Electrodes for Resting EKG & Echo Sound

## FEATURES



### HIGH CHLORIDE WET GEL

High chloride wet gel provides fast and clear traces.



### HELP REDUCE ARTIFACTS

Self-adhering electrodes help reduce muscle artifact or repeat testing.



### AVAILABLE OPTIONS

Available in easy-to-peel vinyl tape or occlusive foam.



### DEFIBRILLATOR COMPATIBLE

Allows for defibrillation if necessary.



### HYPOALLERGENIC

Hypoallergenic adhesive, not made with natural rubber latex.

Vermed SilveRest™ Electrodes set the standard in diagnostic electrodes.

## PACKAGING

### RESEALABLE BULK POUCH

Ziplock pouch contains 30 or more electrodes

### CASE

FDA / GS1 Compliant Labels



### HYGIENIC, DISPOSABLE SINGLE-USE PACKS

Reduces risk of cross contamination of unused electrodes

### INNER BOX

Single use packs are packaged in inner boxes within the case

1.800.669.6905 | [HS.NisshaMedical.com](http://HS.NisshaMedical.com)

Patient Monitoring & Medical Surgical Consumables

Latex-free



Wet Gel



Solid Gel



# SilveRest™ Electrodes

Adult Diagnostic for Resting EKG & Echo Sound by Vermed

Vermed SilveRest™ Electrodes set the industry standard in diagnostic electrodes. Superior Silver/Silver Chloride technology and quality construction produce an electrode that performs every time. Excellent adhesion and clear traces mean no repeat testing.

## A10009 - 1.375" / 34.925mm - ROUND VINYL TAPE ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
A10009	10	10	200	1200		24 Hours
A10009-100	10	100	n/a	1000		

## A10009-100F - 1.375" / 34.925mm - ROUND FOAM ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
A10009-100F	10	100	n/a	1000		24 Hours

## A10013-100S - 1" x 1.5" / 25.4mm x 38.1mm - RECT. FOAM ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel	Duration of Use
A10013-100S	5	100	n/a	1000		24 Hours

## A10049 - 0.8125" x 1.75" / 20.638mm x 44.45mm - RECT. FOAM ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
A10049	5	50	500	2000		24 Hours