



TenderTrode Plus™

Pre-Wired Electrodes & Thermal Probe Covers

FEATURES



GENTLE ON SKIN

Designed specifically for NICU, is gentle and safe for delicate neonatal skin.



SOLID GEL

Leaves no residue on skin and reduces clean-up.



HYDROCOLLOID BORDER

Resists breakdown from moisture or heat.



5 DAY WEAR TIME



Provides long-term attachment with no skin irritation allowing for fewer electrode changes and saving nursing time and resources.

Vermed's TenderTrode Plus™ line of hydrocolloid products are designed for the unique needs of the neonatal market.

PACKAGING

CASE

FDA / GS1 Compliant Labels




HYGIENIC, DISPOSABLE SINGLE-USE PACKS

Reduces risk of cross contamination of unused electrodes

PRE-WIRED ELECTRODE ANATOMY

Each electrode is designed with gentle, hypoallergenic adhesives and soft, conforming materials to ensure comfort and effectiveness on young patients.



**Enlarged for display*

TenderTrode Plus™


Pre-Wired Electrodes & Thermal Probe Covers

Vermed's TenderTrode Plus™ line of hydrocolloid products are designed for the unique needs of the neonatal market. Each product uses a hydrocolloid instead of aggressive tape adhesives to prevent irritation to tender skin while providing long-term attachment. This, coupled with the use of soft conforming materials, ensures patient comfort and product effectiveness in the NICU.



A10055-SRT - 1" / 25.4MM - ROUND PRE-WIRED CLOTH - 3 LEAD - RT

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
A10055-SRT	3	3	n/a	60		 120 Hours



A10077-SRT - 1" / 25.4MM - ROUND PRE-WIRED CLOTH - 5 LEAD - RT

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
A10077-SRT	5	5	n/a	100		 120 Hours



A10058-SRT - 1" X 0.375" / 25.4MM X 9.525MM - RECT. PRE-WIRED CLOTH - RT

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
A10058-SRT	3	3	n/a	60		 120 Hours



A10076 - 1" / 25.4MM - HYDROCOLLOID TEMPERATURE PROBE COVER

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
A10076	1	50	n/a	500		 120 Hours



A10085 - CANNULA HOLDER

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
A10085	1	n/a	n/a	100		 120 Hours



A10086 - MINI CANNULA HOLDER FOR NEONATES

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
A10086	1	n/a	n/a	100		 120 Hours

CLINICAL FEEDBACK

Recent clinical feedback indicates a **98% reduction in faulty incubator and monitor alarms** tripped when standard electrodes peel off in high heat and humidity.

TenderTrode™ Plus is now **leading the way** in neonatal patient care. **Significant cost savings** are evident due to limited replacement, a function of our patented hydrocolloid adhesive border. These electrodes have a 5 day wear time, guaranteed!

Reducing the removal and replacement of standard monitoring electrodes **ensures the body is focused on essential growth and repair**. Reduced monitor alarms ensures the comfort of patients and parents alike while reducing nursing staff workload.